



Wireless connectivity

Bernafon accessories.
Everything that's
needed, perfectly
networked.

With intermediary accessories or direct-to-ear streaming, Bernafon's accessories enable connectivity between hearing aids and from hearing aids to a host of sound sources.





Get control with the Bernafon App

Make the most of your hearing aid experience with the easy-to-use Bernafon App. Control your hearing aids, personalize your sound, and connect to your favorite devices – all discreetly at your fingertips. Get remote counseling from anywhere and set and achieve goals with the Hearing Coach feature. Even use the app to find your hearing aids if they're lost. Now compatible with Apple Watch, for even more discreet control.

Download the FREE app now





SoundClip-A

The all-round device for hands-free phone conversations from iPhone and modern smartphones. Worn around the neck or clipped on the shirt, it makes for improved interpersonal





To check compatibility, see this [list of compatible devices](#).

Compatible with [Alpha XT](#), [Alpha](#), [Viron](#), [Leox](#) and [Zerena](#).



TV-A (TV adapter)

This adapter streams the sound from your TV directly into the hearing aid. Fast and easy to set up, it can pair with an unlimited number of Bernafon hearing aids.

Compatible with [Alpha XT](#), [Alpha](#), [Viron](#), [Leox](#) and [Zerena](#).





RC-A (remote control)

This remote control makes for smooth and discreet hearing aid operation. Lightweight, with excellent tactility and low power consumption.

Compatible with [Alpha XT](#), [Alpha](#), [Viron](#) and [Leox](#) and [Zerena](#).





2.4 GHz Bluetooth[®] wireless technology

Bernafon's hearing aids are Made for iPhone[®] and allow for direct communication and control with an iPhone[®], iPad Pro[®], iPad Air[®], iPad mini[®], and iPod touch[®].

With 2.4 GHz technology, users can enjoy sound streaming directly to their hearing aids. Bernafon Alpha XT and Alpha support ASHA (Audio Streaming for Hearing Aids). With ASHA, you can stream audio from your mobile device directly to your hearing aids.

To check compatibility, see this [list of compatible devices](#).

Compatible with [Alpha XT](#), [Alpha](#), [Viron](#), [Leox](#), and [Zerena](#).





Hands-free communication

Bernafon now offers hands-free communication from iPhone[®] and iPad[®]. Users can leave their phone aside and instead accept, reject, or end a call via hearing aid push button. The hearing aids will pick up their voice and stream sound both ways, leaving their hands free for other tasks.

To check compatibility, see this [list of compatible devices](#).

Compatible with [Alpha XT](#) and [Alpha](#).





SoundGate 3

SoundGate gives you the freedom to listen with your hearing aids to a number of wired and wireless audio sources. It works as an interface between some Bernafon hearing aids and Bluetooth[®] enabled mobile phones.

Compatible with **Saphira** (depending on the wireless capability of the hearing aid style).





SoundGate app

The SoundGate app goes beyond the usual controls of a hearing aid. Use the app on your mobile phone together with the SoundGate to switch programs, turn the volume up or down, and connect to multiple devices.

To check compatibility, see this [list of compatible devices](#).




Compatible with **Saphira** (depending on the wireless capability of the hearing aid style).





Professional

You can also speak to your hearing care professional if you have any questions



[Hearing aid users](#) } [Hearing aids](#) } [Connectivity](#)

[Online hearing test](#)

[Library](#)

[Find a hearing care professional](#)

[Contact](#)

Select country: [United Kingdom](#)

© Copyright Bernafon 2023

[Privacy Notice](#) [Disclaimer](#) [Cookie Policy](#) [Declarations of Conformity](#)

